



# Disposable Heavy Duty Stretcher Sheet with Restraint Belts Access

- Single patient use reduces cross contamination risk
- Made from durable and breathable material for excellent patient comfort
- Fitted corner for secure fitting to standard stretcher mattresses
- Lateral transfer up to 280 kg
- Allows ease of use for most shoulder harness restraints
- Extra length




Instructions:

## Specifications:

|                 |                                |
|-----------------|--------------------------------|
| Product code :  | BL126                          |
| Material :      | 80 GSM Non-Woven Polypropylene |
| Colour :        | White                          |
| Slide Transfer: | up to 280 kg                   |
| Dimensions :    | 80 cm x 200 cm                 |
| Pack Size :     | Pk 50                          |

### A Fitting to the Mattress

Attach the fitted sheet corners to ensure a secure fitting to the mattress.



### B Use as Patient Transfer Sheet

Ensure that the fitting of both ends of the sheet is loosened before attempting lateral transfer of the patient.

Ensure the sheet is not compromised and take a firm grip.



[www.medilinen.com](http://www.medilinen.com)

[info@medilinen.com](mailto:info@medilinen.com)

