



Disposable Heavy Duty Stretcher Sheet


- Single patient use reduces cross contamination risk
- Made from durable and breathable material for excellent patient comfort
- Fluid resistant polypropylene
- Fitted corners for secure fitting to standard stretcher mattresses
- Unsoiled items can be recycled
- Individually packaged



Instructions:

A Fitting to the Mattress

Attach the fitted sheet corners to ensure a secure fitting to the mattress.



B Use as Patient Transfer Sheet

Ensure that the fitting of both ends of the sheet is loosened before attempting lateral transfer of the patient.

Ensure the sheet is not compromised and take a firm grip.



Specifications:

Product code :	BL125-F
Material :	40 gm/m2 Polypropylene
Colour :	White
Slide Transfer Limit:	280kg
Dimensions :	76 cm x 197 cm
Pack Size :	Pk 50

www.medilinen.com

info@medilinen.com

